

Is Your Practice Looking For The Most Efficient BP Monitoring System?

A BpTRU System effectively replaces your 24-hr ambulatory monitor by taking multiple measurements in 6 minutes to eliminate any “white coat” effect on the very first visit. This saves you time and money as well as unnecessary extra visits by patients.

It is also a more accurate alternative to a waiting room monitor, as blood pressure measurements taken in front of other patients are far from ideal.

Save on annual calibration costs and downtime as your BpTRU self-calibrates every time you switch it on guaranteeing 100% reliability.

**BHS A/A
Validated**

Its range of Latex-free cuffs fit arm circumferences up to 52 cm, so you'll always have the right size for accurate measurements.




Secure your BpTRU to the wall and run it from the mains or position it anywhere in the practice on its rolling stand and run it from its rechargeable battery.

Its new software enables you to operate the BpTRU from another room and automatically save readings into your EMR system.



Achieve the ultimate in consistency by replacing the variety of sphygmomanometers in your practice with just one system and immediately streamline all your BP measurements. Asking your patients to report to the BpTRU station 10 minutes before the GP appointment is plenty of time for multiple measurements and data to be saved on file and be on hand ready for the consultation.



Designed by doctors and nurses, for practices - not consumers. It's a professional BP monitor, so you won't find BpTRU on a supermarket shelf. It's also supported by a 3-year guarantee and next-day replacement.



Call 01273 494401 or email bptru@tyrrellhealthcare.com for an info pack, details of purchase and rental options or to arrange a demo of the BpTRU System.

