



Accessory brackets shown above are available (blood pressure cuff holder and temperature probe holder).



Includes basket, C-clamp, pole and base assembly

PRODUCT FEATURES

- ▶ Sturdy aluminum pole
- ▶ 5 dual wheels, 2 with brakes
- ▶ Large basket
- ▶ Includes C-clamp and hardware needed to mount BpTRU™ Vital Sign Monitor
- ▶ Easy-to-assemble

PRODUCT BENEFITS

- ▶ Allows user to transport BpTRU™ unit and all accessories from one location to another with ease and convenience
- ▶ Large basket allows storage for cuffs and accessories
- ▶ Saves time, space and money by protecting the BpTRU™ Vital Signs Monitor from accidents and disappearance
- ▶ Optimal viewing of BpTRU™ Vital Signs Monitor from the standing or seated position

Good blood pressure measurement begins with **quality** accessories

MODEL DESCRIPTION AND ORDERING NUMBERS

BPM-300-45: BpTRU™ Mobile Stand includes a split anodized aluminum column, large black powder coated basket with cord wrap, C-clamp to mount BpTRU™ Vital Signs monitor, sturdy five split wheel caster base assembly (two casters with brakes), hardware and easy-to-assemble instruction sheet.

WARRANTY INFORMATION

VSM MedTech Devices Inc. warrants to the original purchaser that the BpTRU™ will be free from defects in material and workmanship for two years from the date of the original purchase from VSM or its authorized resellers. All BpTRU™ accessories are warranted for a period of one year from the date of the original purchase from VSM or its authorized resellers. Some restrictions apply.

GENERAL INFORMATION

STAND SIZE

Floor to top column height	4' 3" (129.5 cm)
Anodized column diameter	1 1/4" (3.2cm)
Total anodized column length	47" (119.4 cm)
Top column length	23 1/2" (59.7 cm)
Bottom column length	23 1/2" (59.7 cm)

BLACK COMPOSITE BASE ASSEMBLY

Base assembly height	5" (12.7 cm)
Base assembly diameter	23 7/16" (59.6 cm)
Black split wheel caster diameter	2 3/8" (6 cm)

BASKET

Black wrinkle, powder coated basket with cord wrap 12" x 8" x 6" (30.5 x 20.3x 15.2cm)

C-CLAMP

Extruded aluminum clamp fits column diameters of 1" to 2"

WEIGHT

11 lbs with packaging (5 kg)

PACKAGING

24" x 24" x 7" single wall 200lb test corrugated carton (61 x 61 x 17.8 cm)

ACCESSORIES

BPTRU™ NIBP AND TEMPERATURE MONITORS

BPM-100-3	BPM-100-03: BpTRU™ Automated BP Monitor includes power supply, wall-mount bracket and hardware, coiled tubing and three latex-free adult cuffs (small, regular and large).
BPM-200-3	BPM-200-03: BpTRU™ Vital Signs Monitor includes rechargeable battery, USB port, power supply, wall-mount and table-top bracket and hardware, coiled tubing and three latex-free adult cuffs (small, medium, and large).
BPM-300T	BPM-300T: BpTRU™ Vital Signs Monitor includes rechargeable battery, temperature, USB port, power supply, wall-mount and table-top bracket and hardware, coiled tubing and latex-free adult cuffs (small, medium, large). Oral / axillary temperature probe and 2 boxes of probe covers are included for continuous and single (predictive) measurements. Rectal temperature probe is optional.

TEMPERATURE ACCESSORIES AND SUPPLIES

BPM-300-20	Oral/Axillary Temperature Probe (blue) and Adhesive Holders
BPM-300-21	Rectal Temperature Probe (red) and Adhesive Holders
BPM-300-24	Box of Disposable Probe Covers (100 covers, packaged 100/box)
BPM-300-25	Case of Disposable Probe Covers (5000 covers, packaged 100/box)

MOUNTING ACCESSORIES AND SUPPLIES

BPM-300-40	Bracket Assembly for BPM-300 model (wall mount or table-top)
BPM-300-41	Mobile Stand C-clamp
BPM-300-42	Temperature Probe Cover Bracket
BPM-300-43	Blood Pressure Cuff Bracket
BPM-300-45	BpTRU™ MOBILE STAND UNIT (compatible with BPM-200 and BPM-300 models only) includes: <ul style="list-style-type: none"> ▶ Mobile Stand C-clamp and hardware ▶ Large Basket ▶ Pole and Base Assembly