

WatchBP 03

24-Hour ABPM



It's also got the
smallest price:
£795 + VAT

Watch BP 03 Benefits at a glance:

- 24-Hour ABPM with quick-set buttons for 15/30/45/60 minute day/night period measurement intervals for complete convenience.
- 7-Day Mode enables simple collection of blood pressure measurements spread across 7 days in morning and evening for patient convenience, as alternative to multiple readings in 24 hours.
- 250 measurements stored with times and averages.
- BHS A/A Grade and ESH Pass.
- Supplied with medium and large cuffs to fit arm circumferences between 22 and 42 cm.
- Analysis software and USB cable included for quick upload of data. No docking station required
- Lightweight at just 260g including batteries and at just 115 x 80 x 35, is possibly the most comfortable ABPM for patients.
- 3-Year guarantee.
- 4 x 1.5 V AAA batteries supplied.

To order your Watch BP 03, call Microlife's UK Distributor on 01273 494401 or email us on watchbp@tyrrellhealthcare.com www.tyrrellhealthcare.com